

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FOXD4L1 Polyclonal Antibody

多克隆抗体

产品货号: PA04062

第三版

描述: This gene is a member of the forkhead/winged-helix (FOX) family of transcription factors with highly conserved FOX DNA-binding domains. Members of the FOX family of transcription factors are regulators of embryogenesis and may play a role in human cancer. This gene lies in a region of chromosome 2 that surrounds the site where two ancestral chromosomes fused to form human chromosome 2. This region is duplicated elsewhere in the human genome, primarily in subtelomeric and pericentromeric locations, thus multiple copies of this gene have been found. FOXD4L1 (Forkhead Box D4-Like 1) is a Protein Coding gene. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and RNA polymerase II transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is FOXD4.

抗原: Synthesized peptide derived from the C-terminal region of human FoxD4L1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: