

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FOXN4 Polyclonal Antibody

多克隆抗体

产品货号: PA04072

第三版

描述:Members of the winged-helix/forkhead family of transcription factors, such as FOXN4, are characterized by a 110-amino acid DNA-binding domain that can fold into a variant of the helix-turn-helix motif consisting of 3 alpha helices flanked by 2 large loops or wings. These transcription factors are involved in a variety of biologic processes as key regulators in development and metabolism.FOXN4 (Forkhead Box N4) is a Protein Coding gene. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and chromatin binding. An important paralog of this gene is FOXN1.

抗原:Synthesized peptide derived from the C-terminal region of human FoxN4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: