

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human THBS4 Polyclonal Antibody

多克隆抗体

产品货号: PA04075

第三版

描述: THBS4 (Thrombospondin 4) is a Protein Coding gene. Diseases associated with THBS4 include Mycetoma and Myocardial Infarction. Among its related pathways are Notch signaling pathway (KEGG) and Focal Adhesion. GO annotations related to this gene include calcium ion binding and heparin binding. An important paralog of this gene is THBS3. The protein encoded by this gene belongs to the thrombospondin protein family. Thrombospondin family members are adhesive glycoproteins that mediate cell-to-cell and cell-to-matrix interactions. This protein forms a pentamer and can bind to heparin and calcium. It is involved in local signaling in the developing and adult nervous system, and it contributes to spinal sensitization and neuropathic pain states. This gene is activated during the stromal response to invasive breast cancer. It may also play a role in inflammatory responses in Alzheimer's disease. Alternative splicing results in multiple transcript variants.

抗原: Synthesized peptide derived from the Internal region of human Thrombospondin 4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: