

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FOXR1 Polyclonal Antibody

多克隆抗体

产品货号: PA04084

第三版

描述:FOXR1 is a member of the forkhead box (FOX) family of transcription factors. FOX family members are monomeric, helix turn helix proteins with a core DNA binding domain of approximately 110 amino acids. Many FOX transcription factors play roles in determining cell fates during early development. This forkhead box protein lacks the C terminal basic region found in many other FOX family members. It is located within the 11q23.3 region which is commonly deleted in neuroblastomas.

抗原:Synthesized peptide derived from the C-terminal region of human FoxR1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: