

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MMTAG2 Polyclonal Antibody

多克隆抗体

产品货号: PA04090

第三版

描述: C1orf35 (Chromosome 1 Open Reading Frame 35) is a Protein Coding gene. Diseases associated with C1orf35 include Multiple Myeloma. Among its related pathways are Innate Immune System. GO annotations related to this gene include poly(A) RNA binding.

抗原: Synthesized peptide derived from the N-terminal region of human MMTAG2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: