

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FRS2 Polyclonal Antibody

多克隆抗体

产品货号: PA04096

第三版

**描述:**Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.

**抗原:**Synthesized peptide derived from human FRS2 around the non-phosphorylation site of Tyr436.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Highly expressed in heart, brain, spleen, lung, liver, skeletal muscle, kidney and testis.