

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat FSHR Polyclonal Antibody

多克隆抗体

产品货号: PA04098

第三版

**描述:**The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants. Its activation is necessary for the hormonal functioning of FSH. FSHRs are found in the ovary, testis, and uterus.

**抗原:**Synthesized peptide derived from the Internal region of human FSHR

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Sertoli cells and ovarian granulosa cells.