

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat NOLC1 Polyclonal Antibody

多克隆抗体

产品货号: PA04107

---

第三版

**描述:**NOLC1 (Nucleolar And Coiled-Body Phosphoprotein 1) is a Protein Coding gene. Diseases associated with NOLC1 include Treacher Collins Syndrome 1. GO annotations related to this gene include poly(A) RNA binding and GTP binding.

**抗原:**Synthesized peptide derived from the C-terminal region of human Nopp140

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**