

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat HSP60 Polyclonal Antibody

多克隆抗体

产品货号: PA04124

---

第三版

**描述:** HSPs is a group of stress protein cells induced by stress stimulation, and it is related to tumorigenesis, proliferation and differentiation. HSPs can be classified into 3 types according to their different molecular weight, and the distribution and function of each type of HSPs are different. HSP60 is mainly used in research of tongue squamous cell carcinoma, colon cancer, pancreatic cancer and breast cancer, which may have prognostic significance.

**抗原:** Recombinant Protein

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:**