

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Galectin 4 Polyclonal Antibody

多克隆抗体

产品货号: PA04125

第三版

描述:The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The expression of this gene is restricted to small intestine, colon, and rectum, and it is underexpressed in colorectal cancer. LGALS4 (Galectin 4) is a Protein Coding gene. Diseases associated with LGALS4 include Paraneoplastic Cerebellar Degeneration and Colorectal Cancer. Among its related pathways are Cell adhesion_Cell-matrix glycoconjugates. GO annotations related to this gene include carbohydrate binding. An important paralog of this gene is LGALS9.

抗原:Synthesized peptide derived from the Internal region of human Galectin-4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: