

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Galectin 9 Polyclonal Antibody

多克隆抗体

产品货号: PA04127

第三版

描述:The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The protein encoded by this gene is an S-type lectin. It is overexpressed in Hodgkin's disease tissue and might participate in the interaction between the H&R;S; cells with their surrounding cells and might thus play a role in the pathogenesis of this disease and/or its associated immunodeficiency. Multiple alternatively spliced transcript variants have been found for this gene.LGALS9 (Galectin 9) is a Protein Coding gene. Among its related pathways are RET signaling and Innate Immune System. GO annotations related to this gene include signal transducer activity and galactose binding. An important paralog of this gene is LGALS9B.

抗原:Synthesized peptide derived from the Internal region of human Galectin-9

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Peripheral blood leukocytes and lymphatic tissues. Overexpressed in Hodgkin disease tissue.