

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GAS3 Polyclonal Antibody

多克隆抗体

产品货号: PA04131

第三版

描述: This gene encodes an integral membrane protein that is a major component of myelin in the peripheral nervous system. Studies suggest two alternately used promoters drive tissue-specific expression. Various mutations of this gene are causes of Charcot-Marie-Tooth disease Type 1A, Dejerine-Sottas syndrome, and hereditary neuropathy with liability to pressure palsies. Alternative splicing results in multiple transcript variants. PMP22 (Peripheral Myelin Protein 22) is a Protein Coding gene. Diseases associated with PMP22 include Neuropathy, Recurrent, With Pressure Palsies and Charcot-Marie-Tooth Disease, Type 1A. Among its related pathways are Neural Crest Differentiation and $\alpha 6 \beta 1$ and $\alpha 6 \beta 4$ Integrin signaling. An important paralog of this gene is EMP2.

抗原: Synthesized peptide derived from the C-terminal region of human GAS3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20°C 保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 $\mu\text{g/ml}$.

特异性: