

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GSDMC Polyclonal Antibody

多克隆抗体

产品货号: PA04135

第三版

描述:GSDMC (Gasdermin C) is a Protein Coding gene. An important paralog of this gene is GSDMA.

抗原:Synthesized peptide derived from the Internal region of human Gasdermin C

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Expressed mainly in trachea and spleen (PubMed:11223543). In the esophagus, expressed in differentiating cells and probably in differentiated cells. Also detected in gastric epithelium (PubMed:19051310).