
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat EPG5 Polyclonal Antibody

多克隆抗体

产品货号: PA04141

第三版

描述: This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions. Mutations in this gene cause Vici syndrome. EPG5 (Ectopic P-Granules Autophagy Protein 5 Homolog) is a Protein Coding gene. Diseases associated with EPG5 include Vici Syndrome and Spastic Paraplegia 49, Autosomal Recessive.

抗原: Recombinant Protein

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: