

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human UBAC1 Polyclonal Antibody

多克隆抗体

产品货号: PA04145

第三版

描述:KPC2 (Kip1 ubiquitylation-promoting complex 2) together with KPC1 forms the ubiquitin ligase KPC, which regulates degradation of the cyclin-dependent kinase inhibitor p27 at the G1 phase of the cell cycle. KPC2 contains a ubiquitin-like (UBL) domain, two ubiquitin-associated (UBA) domains, and a heat shock chaperonin-binding (ST11) domain.

抗原:Synthesized peptide derived from the Internal region of human GBDR1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.