

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GBP3 Polyclonal Antibody

多克隆抗体

产品货号: PA04148

第三版

描述: This gene encodes a member of the guanylate-binding protein (GBP) family. GBPs specifically bind guanine nucleotides (GMP, GDP, and GTP) and contain two of the three consensus motifs found in typical GTP-binding proteins. The encoded protein interacts with a member of the germinal center kinase family. Multiple transcript variants encoding different isoforms have been found for this gene. GBP3 (Guanylate Binding Protein 3) is a Protein Coding gene. Among its related pathways are Interferon gamma signaling and Innate Immune System. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is GBP1.

抗原: Synthesized peptide derived from the C-terminal region of human GBP3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: