

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GBP4 Polyclonal Antibody

多克隆抗体

产品货号: PA04150

第三版

描述:Guanylate-binding proteins, such as GBP4, are induced by interferon and hydrolyze GTP to both GDP and GMP.GBP4 (Guanylate Binding Protein 4) is a Protein Coding gene. Among its related pathways are Interferon gamma signaling and Innate Immune System. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is GBP7.

抗原:Synthesized peptide derived from the C-terminal region of human GBP4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: