

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GCP3 Polyclonal Antibody

多克隆抗体

产品货号: PA04159

第三版

描述:TUBGCP3 (Tubulin Gamma Complex Associated Protein 3) is a Protein Coding gene. Diseases associated with TUBGCP3 include Cardiomyopathy, Dilated, 1P. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Chks in Checkpoint Regulation. GO annotations related to this gene include structural molecule activity and gamma-tubulin binding. An important paralog of this gene is TUBGCP2.

抗原:Synthesized peptide derived from the Internal region of human GCP3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed.