

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GCP5 Polyclonal Antibody

多克隆抗体

产品货号: PA04163

第三版

**描述:**TUBGCP5 (Tubulin Gamma Complex Associated Protein 5) is a Protein Coding gene. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Chks in Checkpoint Regulation. GO annotations related to this gene include microtubule binding. An important paralog of this gene is TUBGCP6.

**抗原:**Synthesized peptide derived from the Internal region of human GCP5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Widely expressed, with highest levels in heart and skeletal muscle and moderate levels in brain.