

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GCP6 Polyclonal Antibody

多克隆抗体

产品货号: PA04164

第三版

描述:The protein encoded by this gene is part of a large multisubunit complex required for microtubule nucleation at the centrosome. TUBGCP6 (Tubulin Gamma Complex Associated Protein 6) is a Protein Coding gene. Diseases associated with TUBGCP6 include Microcephaly And Chorioretinopathy, Autosomal Recessive, 1 and Autosomal Recessive Chorioretinopathy-Microcephaly Syndrome. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Chks in Checkpoint Regulation. GO annotations related to this gene include microtubule binding. An important paralog of this gene is TUBGCP2.

抗原:Synthesized peptide derived from the Internal region of human GCP6

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: