

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NMBR Polyclonal Antibody

多克隆抗体

产品货号: PA04165

第三版

描述: This gene encodes a 7-transmembrane G protein-coupled receptor that binds neuromedin B, which is a growth factor and mitogen for gastrointestinal epithelial tissue and for normal and neoplastic lung. This receptor may play a role in smooth muscle contraction, neuronal responses, and the regulation of cell growth. Antagonists of this receptor have a potential therapeutic use in inhibiting tumor cell growth. Polymorphisms in this gene may be associated with a susceptibility for schizophrenia. Alternative splicing of this gene results in multiple transcript variants. NMBR (Neuromedin B Receptor) is a Protein Coding gene. Diseases associated with NMBR include Lung Cancer and Small Cell Cancer Of The Lung, Somatic. Among its related pathways are Peptide ligand-binding receptors and RET signaling. GO annotations related to this gene include G-protein coupled receptor activity and bombesin receptor activity. An important paralog of this gene is GRPR.

抗原: Synthetic Peptide

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in epididymis (at protein level).