

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MMP9 Polyclonal Antibody

多克隆抗体

产品货号: PA04177

第三版

描述: May play an essential role in local proteolysis of the extracellular matrix and in leukocyte migration. Could play a role in bone osteoclastic resorption. Cleaves KiSS1 at a Gly-Leu bond. Cleaves type IV and type V collagen into large C-terminal three quarter fragments and shorter N-terminal one quarter fragments. Degrades fibronectin but not laminin or Pz-peptide.

抗原: Synthesized peptide derived from the C-terminal region of human MMP-9

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Produced by normal alveolar macrophages and granulocytes.