

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GFAP Polyclonal Antibody

多克隆抗体

产品货号: PA04182

第三版

描述: This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

抗原: Synthesized peptide derived from human GFAP around the non-phosphorylation site of Ser38.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in cells lacking fibronectin.