

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GIT1 Polyclonal Antibody

多克隆抗体

产品货号: PA04194

第三版

**描述:** GIT1 (GIT ArfGAP 1) is a Protein Coding gene. Diseases associated with GIT1 include Chromosome 17Q11.2 Deletion Syndrome, 1.4-Mb and Chromosome 17Q11.2 Deletion Syndrome. Among its related pathways are Arf6 signaling events and Adhesion. GO annotations related to this gene include GTPase activator activity. An important paralog of this gene is GIT2.

**抗原:** Synthesized peptide derived from the C-terminal region of human GIT1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**