

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GK2 Polyclonal Antibody

多克隆抗体

产品货号: PA04203

第三版

**描述:**GK2 (Glycerol Kinase 2) is a Protein Coding gene. Diseases associated with GK2 include Glycerol Kinase Deficiency. Among its related pathways are triacylglycerol biosynthesis and Metabolism. GO annotations related to this gene include phosphotransferase activity, alcohol group as acceptor and glycerol kinase activity. An important paralog of this gene is GK.

**抗原:**Synthesized peptide derived from the Internal region of human GK2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**