

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR100 Polyclonal Antibody

多克隆抗体

产品货号: PA04238

第三版

描述:High affinity receptor for INSL5. Also acts as receptor for RLN3/relaxin-3, as well as bradykinin and kallidin. Binding of the ligand inhibit cAMP accumulation.

抗原:Synthesized peptide derived from the C-terminal region of human GPR100

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in a broader range of tissues including brain, kidney, testis, thymus, placenta, prostate, salivary gland, thyroid and colon.