
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR119 Polyclonal Antibody

多克隆抗体

产品货号: PA04257

第三版

描述: Receptor for the endogenous fatty-acid ethanolamide oleoylethanolamide (OEA) and lysophosphatidylcholine (LPC). Functions as a glucose-dependent insulinotropic receptor. The activity of this receptor is mediated by G proteins which activate adenylate cyclase. Seems to act through a G(s) mediated pathway.

抗原: Synthesized peptide derived from the Internal region of human GPR119

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Predominantly expressed in the pancreas, especially in the islets.