

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DREG Polyclonal Antibody

多克隆抗体

产品货号: PA04261

第三版

**描述:**Orphan receptor. May be required for normal differentiation of promyelinating Schwann cells and for normal myelination of axons. Signals probably through G-proteins to transiently elevate cAMP levels.

**抗原:**Synthesized peptide derived from the Internal region of human GPR126

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in placenta and to a lower extent in pancreas and liver. Detected in aortic endothelial cells but not in skin microvascular endothelial cells.