

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CABP2 Polyclonal Antibody

多克隆抗体

产品货号: PA04263

第三版

描述:CABP2 (Calcium Binding Protein 2) is a Protein Coding gene. Diseases associated with CABP2 include Deafness, Autosomal Recessive 93 and Dfnb93 Nonsyndromic Hearing Loss And Deafness. Among its related pathways are Visual Cycle in Retinal Rods. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CABP1. This gene belongs to a subfamily of calcium binding proteins that share similarity to calmodulin. Like calmodulin, these family members can likely stimulate calmodulin-dependent kinase II and the protein phosphatase calcineurin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Two transcript variants encoding different isoforms have been found for this gene.

抗原:Synthetic Peptide

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Retina.