

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR172A Polyclonal Antibody

多克隆抗体

产品货号: PA04277

第三版

描述:Riboflavin transporter. Riboflavin transport is Na(+)-independent but moderately pH-sensitive. Activity is strongly inhibited by riboflavin analogs, such as lumiflavin. Weakly inhibited by flavin adenine dinucleotide (FAD) and flavin mononucleotide (FMN). In case of infection by retroviruses, acts as a cell receptor to retroviral envelopes similar to the porcine endogenous retrovirus (PERV-A).

抗原:Synthesized peptide derived from the Internal region of human GPR172A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in brain, fetal brain and salivary gland. Weakly expressed in other tissues.