

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GPR180 Polyclonal Antibody

多克隆抗体

产品货号: PA04289

第三版

**描述:** This gene encodes a protein that is a member of the G protein-coupled receptor superfamily. This protein is produced predominantly in vascular smooth muscle cells and may play an important role in the regulation of vascular remodeling. GPR180 (G Protein-Coupled Receptor 180) is a Protein Coding gene. Diseases associated with GPR180 include La Crosse Encephalitis and Fungal Meningitis.

**抗原:** Synthesized peptide derived from the Internal region of human GPR180

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**