

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR62 Polyclonal Antibody

多克隆抗体

产品货号: PA04308

第三版

描述:GPR62 (G Protein-Coupled Receptor 62) is a Protein Coding gene. Among its related pathways are GPCRs, Other. GO annotations related to this gene include G-protein coupled receptor activity. An important paralog of this gene is GPR61.

抗原:Synthesized peptide derived from the Internal region of human GPR62

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in brain; detected in the basal forebrain, frontal cortex, caudate, putamen, thalamus and hippocampus.