

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DREAM Polyclonal Antibody

多克隆抗体

产品货号: PA04314

第三版

描述: This gene encodes a member of the family of voltage-gated potassium (Kv) channel-interacting proteins, which belong to the recoverin branch of the EF-hand superfamily. Members of this family are small calcium binding proteins containing EF-hand-like domains. They are integral subunit components of native Kv4 channel complexes that may regulate A-type currents, and hence neuronal excitability, in response to changes in intracellular calcium. The encoded protein also functions as a calcium-regulated transcriptional repressor, and interacts with presenilins. Alternatively spliced transcript variants encoding different isoforms have been described. KCNIP3 (Potassium Voltage-Gated Channel Interacting Protein 3) is a Protein Coding gene. Diseases associated with KCNIP3 include Alzheimer Disease. Among its related pathways are Cardiac conduction and Development Ligand-independent activation of ESR1 and ESR2. GO annotations related to this gene include calcium ion binding and transcription corepressor activity. An important paralog of this gene is KCNIP4.

抗原: Synthetic Peptide

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in brain. Widely expressed at lower levels. Expression levels are elevated in brain cortex regions affected by Alzheimer disease.