

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPRIN1 Polyclonal Antibody

多克隆抗体

产品货号: PA04333

第三版

描述:GPRIN1 (G Protein Regulated Inducer Of Neurite Outgrowth 1) is a Protein Coding gene. Diseases associated with GPRIN1 include Pustulosis Of Palm And Sole and Psoriasis 14, Pustular. Among its related pathways are MECP2 and Associated Rett Syndrome. GO annotations related to this gene include phosphoprotein binding. An important paralog of this gene is GPRIN3.

抗原:Synthesized peptide derived from the Internal region of human GRIN1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed in the central nervous system, with highest levels in spinal cord.