
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GRK1 Polyclonal Antibody

多克隆抗体

产品货号: PA04337

第三版

描述:Phosphorylates rhodopsin thereby initiating its deactivation. This rapid desensitization is essential for scotopic vision and permits rapid adaptation to changes in illumination.

抗原:Synthesized peptide derived from human GRK 1 around the non-phosphorylation site of S21.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Retinal-specific. Expressed in rods and cones cells.