

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FAT10 Polyclonal Antibody

多克隆抗体

产品货号: PA04353

第三版

描述:UBD (Ubiquitin D) is a Protein Coding gene. Diseases associated with UBD include Hepatocellular Carcinoma. Among its related pathways are TRAF Pathway and fMLP Pathway. GO annotations related to this gene include proteasome binding. An important paralog of this gene is UBB.

抗原:Synthesized peptide derived from the Internal region of human FAT10.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Constitutively expressed in mature dendritic cells and B cells. Mostly expressed in the reticuloendothelial system (e.g. thymus, spleen), the gastrointestinal system, kidney, lung and prostate gland.