
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANKRD20A1 Polyclonal Antibody

多克隆抗体

产品货号: PA04355

第三版

描述:ANKRD20A1 (Ankyrin Repeat Domain 20 Family Member A1) is a Protein Coding gene. Diseases associated with ANKRD20A1 include Behcet Syndrome. An important paralog of this gene is ANKRD20A2.

抗原:Synthesized peptide derived from the Internal region of human A20.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: