

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NCAPH Polyclonal Antibody

多克隆抗体

产品货号: PA04363

第三版

描述: This gene encodes a member of the barr gene family and a regulatory subunit of the condensin complex. This complex is required for the conversion of interphase chromatin into condensed chromosomes. The protein encoded by this gene is associated with mitotic chromosomes, except during the early phase of chromosome condensation. During interphase, the protein has a distinct punctate nucleolar localization. Alternatively spliced transcript variants encoding different proteins have been described.

抗原: Synthesized peptide derived from the Internal region of human hCAP-H

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed at low level. Expressed in proliferating cells.