
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CCNB1IP1 Polyclonal Antibody

多克隆抗体

产品货号: PA04376

第三版

描述: E3 ubiquitin-protein ligase. Modulates cyclin B levels and participates in the regulation of cell cycle progression through the G2 phase. Overexpression causes delayed entry into mitosis.

抗原: Synthesized peptide derived from the C-terminal region of human HEI10

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in heart. Detected at intermediate levels in liver and kidney, and at low levels in placenta, brain and lung.