

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HES6 Polyclonal Antibody

多克隆抗体

产品货号: PA04381

第三版

描述: HES6 is a transcription cofactor that itself does not bind to DNA but suppresses both HES1-mediated N box-dependent transcriptional repression and binding of HES1 to E box sequences. It also suppresses HES1-mediated inhibition of the heterodimer formed by ASCL1/MASH1 and TCFEA2/E47, allowing ASCL1 and TCFEA2 to upregulate transcription in its presence. HES6 appears to have novel functions in tumour and cancer cell lines, although these have not yet been elucidated fully.

抗原: Synthesized peptide derived from the N-terminal region of human HES-6

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: