
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Histone 1.0 Polyclonal Antibody

多克隆抗体

产品货号: PA04391

第三版

描述: Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.

抗原: Synthesized peptide derived from the Internal region of human Histone 1.0

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: