

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human VTN Polyclonal Antibody

多克隆抗体

产品货号: PA04415

第三版

描述:The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

抗原:Synthesized peptide derived from the Internal region of human Vitronectin.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Plasma.