

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TDGF1P3 Polyclonal Antibody

多克隆抗体

产品货号: PA04429

第三版

描述:TDGF1P3 (Teratocarcinoma-Derived Growth Factor 1 Pseudogene 3) is a Pseudogene. Diseases associated with TDGF1P3 include Teratocarcinoma. Among its related pathways are Signaling by NODAL and Developmental Biology. GO annotations related to this gene include growth factor activity. The TDGF3 locus has characteristics of a retrotransposon, including lack of introns and a poly(A) sequence.

抗原:Synthesized peptide derived from the Internal region of human Cripto-3.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed weakly in lung, colon and breast. Expressed also strongly in primary cancer tissues; lung and colon cancers.