
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HLX Polyclonal Antibody

多克隆抗体

产品货号: PA04434

第三版

描述:HLX (H2.0 Like Homeobox) is a Protein Coding gene. Diseases associated with HLX include Ovarian Primitive Germ Cell Tumor and Ovarian Endodermal Sinus Tumor. Among its related pathways are IL12-mediated signaling events. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is DBX1.

抗原:Synthesized peptide derived from the C-terminal region of human HLX1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Low level in normal B and T-cells, high level in activated lymphocytes and monocytes. Also found in thymus, tonsil, bone marrow, developing vessels, and fetal brain.