

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Acetyl-TAL1 (Lys221/Lys222/Lys36/Lys37) Polyclonal Antibody

多克隆抗体

产品货号: PA04440

第三版

描述:TAL1 (TAL BHLH Transcription Factor 1, Erythroid Differentiation Factor) is a Protein Coding gene. Diseases associated with TAL1 include Leukemia, Acute Lymphoblastic and T-Cell Acute Lymphocytic Leukemia-1. Among its related pathways are NF-kappaB Signaling and MAPK-Erk Pathway. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and chromatin binding. An important paralog of this gene is LYL1.

抗原:Synthesized peptide derived from human TAL1/2 around the acetylation site of K221.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Leukemic stem cell.