

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HORMAD1 Polyclonal Antibody

多克隆抗体

产品货号: PA04452

第三版

**描述:**This gene encodes a HORMA domain-containing protein. HORMA domains are involved in chromatin binding and play a role in cell cycle regulation. The encoded protein may play a role in meiosis, and expression of this gene is a potential marker for cancer. A pseudogene of this gene is located on the long arm of chromosome 6. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. HORMAD1 (HORMA Domain Containing 1) is a Protein Coding gene. Diseases associated with HORMAD1 include Splenic Infarction. An important paralog of this gene is HORMAD2.

**抗原:**Synthesized peptide derived from the C-terminal region of human HORMAD1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Testis-specific. Over-expressed in carcinomas.