

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HOXA6 Polyclonal Antibody

多克隆抗体

产品货号: PA04455

第三版

描述:In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which Human, Mouse regulate gene expression, morphogenesis, and differentiation.

抗原:Synthesized peptide derived from the Internal region of human HoxA6

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: