

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CALB2 Polyclonal Antibody

多克隆抗体

产品货号: PA04460

第三版

描述:CALB2 (Calbindin 2) is a Protein Coding gene. Diseases associated with CALB2 include Adenofibroma and Adenomatoid Tumor. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CALB1. This gene encodes an intracellular calcium-binding protein belonging to the troponin C superfamily. Members of this protein family have six EF-hand domains which bind calcium. This protein plays a role in diverse cellular functions, including message targeting and intracellular calcium buffering. It also functions as a modulator of neuronal excitability, and is a diagnostic marker for some human diseases, including Hirschsprung disease and some cancers. Alternative splicing results in multiple transcript variants.

抗原:Synthesized peptide derived from the Internal region of human CALB2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Brain.