

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Acetyl-E2F3 (Lys168) Polyclonal Antibody

多克隆抗体

产品货号: PA04461

第三版

**描述:**E2F3 (E2F Transcription Factor 3) is a Protein Coding gene. Diseases associated with E2F3 include Retinoblastoma and Sensory System Cancer. Among its related pathways are Validated targets of C-MYC transcriptional activation and Cellular Senescence. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and core promoter binding. An important paralog of this gene is E2F2. This gene encodes a member of a small family of transcription factors that function through binding of DP interaction partner proteins. The encoded protein recognizes a specific sequence motif in DNA and interacts directly with the retinoblastoma protein (pRB) to regulate the expression of genes involved in the cell cycle. Altered copy number and activity of this gene have been observed in a number of human cancers. There are pseudogenes for this gene on chromosomes 2 and 17. Alternative splicing results in multiple transcript variants.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**