

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TSC1 Polyclonal Antibody

多克隆抗体

产品货号: PA04483

第三版

描述:TSC1 (Tuberous Sclerosis 1) is a Protein Coding gene. Diseases associated with TSC1 include Tuberous Sclerosis-1 and Lymphangioliomyomatosis. Among its related pathways are Development IGF-1 receptor signaling and RET signaling. GO annotations related to this gene include binding and chaperone binding. This gene encodes a growth inhibitory protein thought to play a role in the stabilization of tuberin. Mutations in this gene have been associated with tuberous sclerosis. Alternative splicing results in multiple transcript variants.

抗原:Synthesized peptide derived from Hamartin

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in skeletal muscle, followed by heart, brain, placenta, pancreas, lung, liver and kidney. Also expressed in embryonic kidney cells.